



Equivizor.co.uk



Important information

We hope you are happy with your new **Recovery mask**, we would like to pass on some tips which may help you and allow you to take advantage of all the features and benefits this mask provides.

We believe our product is one of the best on the current market although we know our mask will not stop every horse Houdini out there we have sure helped a lot of happy owners put their Houdini's out of the Fly mask disappearing act business. We are always open to feedback from our customers so we can continue to improve our service to all our clients.

The **Recovery Equivizor** is for horses suffering Uveitis, recovering from eye surgery or habitual eye scratching. Made of soft, Light-UV-tinted PVC. The mask is available in sizes foal, pony, cob, full and extra-large, dark or light tinted, and includes a UV-resistant fly mesh for extra protection if required.

The Recovery mask is made from clear PVC the unique design allows sufficient space in front and around the eyes. Its 3 strategically located straps allows for a secure fit along with the soft padding surrounding the mask for comfort. The EQUIVIZOR Recovery has a mesh mask which can be used on top to protect against insects and provides efficient shading.

This mask is rather handy as for example if a vet needs to look at or medicate your horses eye they can do so by just releasing the 2 lower strap clips and lifting the mask, then just place it back down and do up the clips again. Another great feature of this mask is that it has great ventilation so it's not a warm, dark, moist environment that could lead to bacteria growth with an eye condition or after surgery.

The EQUIVIZOR Recovery is available in a tinted version which blocks 92% of the UV rays which is the maximum blocking factor obtainable. The shading factor is 32%. In combination with the mesh it ensures maximum protection for the eyes.



Now you have learnt a little more about your product let me tell you some rather important tips to enable easy fitting and care of your Equivizor Fly mask.

- Please don't soak the vizor in any solution for an extended amount of time, you can wash it with mild liquid soap and warm water, rinse and wipe dry with a soft cloth. (Please make sure to rinse the black Sherpa lining very well and not leave soapy residue which could get in your horse's eye)
- Due to its secure fit, lack of restrictions and interference to the horse you can ride in the Equivizor Recovery Mask.
- The mask will inevitably scratch over time. How much will depend on the activity of your horse and their environment. It's a lot better the vizor takes the bumps rather than your horse.
- Its advised that you trim the excess straps off and seal them so they do not fray. If you do not trim them they then become great buddy toys for another horse to play with and pull at.
- If you plan on riding your horse in the vizor then leave the straps a little longer so you can lengthen them when wearing a bridle.
- Make sure you do routine checks of your horses' eyes and face to catch any irritation areas sooner rather than later
- The clips are created to be quick release, if your horse was to get caught on something the clips will break away if your horse pulls hard (Replacement clips are available)
- When measuring your horse for fitting please be sure to visit our website and follow our fitting and sizing guide which can be found on our website
- We are usually available via email if you wish to contact us about fitting and cleaning.

^^ Please check fences and the areas your horse is in, check fence line on a regular basis for sharp edges. A little tip if you find sharp edges that your horse could possibly get caught on is to put black duct tape over the damaged section. Also, be sure to check your feeders and water drinkers for any sharp points or protruding areas ^^

In an effort to provide excellent customer service, we would welcome you to send us pictures of your mask in use. this could help us help you with the fitting of your mask and allow us to hear your thoughts on our product.

(please send at least 1 front view and one each side, while taking the photos please make sure they are as parallel as possible so we can see it at the correct angle)

Please find us on Facebook Equivizor UK, if you are satisfied with our product we always appreciate a review on our social media account and on our website.

Please don't hesitate to contact us with any questions or concerns

Email: sales@equivizor.co.uk

Phone : 01323 831710

www.equivizor.co.uk